

Gislars Canyon JOCKO-ET

CH 120.148.466.260 Born: 18.11.2019



Breed: Brown Swiss
Breeder: -

CANYON DE 08 16637282 133/115/+1083-0,39-0,11	CADURA DE 08 16074070	CADENCE	gGZW 110
	LIESE DE 08 15943923	GS HUVI	
JANNA CH 120.123.170.793 94/96/-424+0,12+0,16 3/2 8118-3,97-3,98-646 3. 9419-3,74-3,90-719	ARROW IT 021002045985	ANIBAL	MW 102 FW 96 FIT 112 Conf. 111
	JUDY CH 120.076.876.278 3/3 9837-4,19-3,60-767	NESTA	

BREEDING VALUES (03.04.2024) GYB gGZW 110 (76)

GZW -1, MW -1, FW +0, FIT -1, EXT +0

MILK -82 -0,01 -4 +0,13 +7 **MW 102 (84)**

BEEF **FW 96 (63)**

Net daily gain: 97 (64) **Carcass percentage:** 102 (60) **EUROP trade class:** 95 (60)

FITNESS **ÖZW 114 (80)** **FIT 112 (77)**

Longevity: 108 (67)	Udder health index: 103 (79)	Somatic cell count: 99 (76)	Mastitis: 109 (52)
Persistency: 111 (74)	Fertility index: 112 (65)	Claw health index: 95 (63)	Early fert. disorders: 107 (59)
Yield improvement: 111 (71)	Calving ease direct: 98 (83)	Calving ease maternal: 92 (68)	Cysts: 109 (59)
Milking speed: 106 (80)	Milking behaviour: 79 (53)	Calf vitality: 103 (66)	Milk fever: -

CONFORMATION **0 daughters: 99 - 105 - 110 - 115 - Conf. 111 (76)**

Merkmal	76	88	100	112	124	136
Frame 99						
Rump 105						
Feet & Legs 110						
Udder 115						
Final score 111						
Muscularity 93			weak			heavy
Cross height 103			small			large
Chest width 96			narrow			wide
Body depth 98			shallow			deep
Backline 97			weak			strong
Rump length 106			short			long
Rump width 102			narrow			wide
Rump angle 101			ascending			sloped
Thurl position 100			in the back			in the centre
Rear legs side view 101			straight			sickled
Hock quality 108			swollen			dry
Pasterns 110			weak			strong
Hoof height 107			low			high
Fore udder length 104			short			long
Rear udder width 94			narrow			wide
Rear udder height 109			low			high
Central ligament 101			weak			strong
Udder depth 117			deep			high
Fore udder attachment 103			loose			tight
Udder balance 117			staged			inclined
Teat length 91			short			long
Teat thickness 92			thin			thick
Front teat placement 114			wide			close
Rear teat placement 108			wide			close
Teat direction 115			outwards			inwards
Udder cleanness 101			add. teats			clean

Free of: ARF, SDF, SMF, WEF

